

UNIT-II Instrumental Method of Analysis

In-SEM Exam = 15 Marks

* Question Paper Pattern - (No. Numericals)

All questions are based on Unit-II :-

Ques. 3 (a) — 5 Marks OR Ques. 4 (a) — 5 Marks

(b) — 4 Marks (b) — 4 Marks

(c) — 3 Marks (c) — 3 Marks

(d) — 3 Marks (d) — 3 Marks

Only Theory Questions.

* Reference books for Unit-II :-

(i) Analytical Chemistry - S.M. Khopkar *

(ii) Engg. Chemistry - O.G. Palana *

(iii) Engg. Chemistry - Shashi Chawla

• problem solving •

• theory & practical •

• MCQs & Definitions •

• Revision Notes •

• Air (iv) 830 or H2 (vi) 92 or A = 60

Course Objective -

To acquire the knowledge of electro-analytical techniques that facilitates rapid & precise understanding of material

• To provide the fundamental concepts of analytical chemistry

• To understand the basic principles of various analytical methods

Index / Mind Map of Unit-II

Ch 1: Prerequisite of Electrochemistry -

- Electrochemical & Electrolytic cell.

- Electrode Potential, Nernst Eqn.

(d)

(d) (d) (d)

* Reference Electrode :-

Definition :- The electrode which has stable & reproducible electrode potential & acts as half cell, thereby completing the cell is known as 'reference electrode'.

Types :- Ref. Electrode


Primary
eg. SHE

Secondary

eg. Calomel, Ag-Ag,
Quinn hydrone,

Calomel Electrode :-

Diagram :-




Principle :- It works as an reference electrode whose electrode potential is stable & reproducible. It acts a half cell & consists of Hg, Hg₂Cl₂, KCl, etc.


Construction of a calomel electrode :-

- It consists of a thick glass tube.
- At the bottom of the tube pure Hg is placed, above which paste of $Hg + Hg_2Cl_2$ is placed.
- In which a Pt wire is fitted in a thin glass tube, & the wire is connected to electrometer.
- Remaining glass tube is filled with saturated KCl soln.
- The electrode is connected to salt bridge & whole assemble is sealed.

Representation of Calomel Electrode -



Reaction of calomel electrode -



Equation of calomel electrode -

$$E_{cal} = E_{cal}^\circ - 0.0591 \log [Cl^-]$$

Electrode Potential -

Conc. of KCl	1 N	0.1 N	Saturated
Electrode Pot.	0.2810 V	0.3334 V	0.2422 V

Disadvantages :- It does not work in the soln having temp. more than $50-80^\circ C$ because paste of $Hg + Hg_2Cl_2$ decomposes in pure Hg, thereby disturbing the whole assembly.

Advantages :- Portable, easy to handle, simple equipment, compact.

* Ion selective electrodes (ISE) :-

Definition / Principle :- The electrode which responds to specific ions of the solution, thereby showing particular potential is called as 'ion selective electrodes' or 'chemical sensors'.

Types of ISE →


- i) Glass membrane electrode / Glass electrode
- ii) Solid state memb. electrode.
- iii) Liquid - Liquid membrane electrode.
- iv) Enzyme based membrane / Biochemical electrode
- v) Gas sensing electrodes.

i) Glass Membrane / Glass electrode :- (Indicator Electrode) →


Principle :- When two solutions of different H^+ ions are separated by thin glass membrane a potential difference is developed on the surface of membrane.

Glass membrane acts as an ion exchange i.e. exchange of Na^+ from glass with H^+ from solutions.

Diagram :-



Representation :-



$$\text{Equation: } E_g = E_g^{\circ} + 0.059 / \text{pH}$$

Construction :-

- It consists of a thick glass tube.
- At the bottom of glass tube, there is a thin glass bulb (thin walled).
- ~~A Pt or Ag~~ ^{coated} wire.
- A Pt wire or Ag coated Ag wire is dipped in 0.1 M HCl soln filled in the bulb.
- The whole assembly is sealed.

Merits :-

- i) Portable
- ii) Easy to handle
- iii) Simple equipment.

ii) Solid state membrane electrode :-
This contains homogeneous mixture of a polycrystalline ionic compounds.

Two types - (a) Fluoride ion determination.
(b) Chloride ion determination.


(a) Fluoride ion determination :-

~~Construction / Composition :-~~ The membrane is made from lanthanum tri-fluoride (LaF_3) crystal doped with europium fluoride (EuF_2).

The crystal is sealed at bottom of polymeric container containing internal ref. soln ($\text{NaF} + \text{NaCl}$ or $\text{KF} + \text{KCl}$) with a reference electrode.

Working :- EuF_2 produces holes in the crystal lattice of LaF_3 through which F^- ions can pass. When electrode is in contact with the sample soln, a potential develops across the membrane.

Diagram :-




(b) Chloride ion determination :-

Composition / Construction :- Pallet of Ag_2S & AgCl is used in this electrode using epoxy resin, electrode is made up of teflon.

Working :- When the electrode membrane is in contact with a solution containing chloride ion, an electrode potential develops. This potential is measured against a const. ref. potential.

Diagram :-




iii) Enzyme based membrane electrode :-

Construction / Composition :-

These electrodes use enzymes to convert substances in the solution into ionic products. A mobilized enzyme is applied at the surface of electrode.

Diagram :-




Construction / Composition :-

The enzyme urease is incorporated into a polyacrylamide gel which is allowed to set on the bulb of a glass electrode & filled by a nylon gauze.

Working :-

When the enzyme based membrane elect. is immersed into a solution containing urea NH_4^+ ions are produced which diffuse through the gel, thereby developing a boundary potential.



iv) Gas Sensing Membrane Electrodes :-

These electrodes are specific towards the ions into which gas is converted after it passes through the outer membrane.

e.g. NH_3 , NO_2 , SO_2 , CO_2 & H_2S .

Construction / Composition -

The gas in the test solution diffuses through the membrane & reacts with the internal filling solution to form the ions. These ions are detected using a gas sensing electrode.

Membrane used - hydrophobic gas-permeable porous polymer membrane (silicon rubber/teflon)

Working :- e.g.

e.g. For NH_3 gas determination -

When electrode is dipped in dissolved NH_3 sample soln, it diffuses through memb. thereby forming boundary potential.

The boundary potential is developed betⁿ ammonium ions formed inside & outside ion selective electrode.





Diagram:



* Conductometric :-

- Important terms involved in conductometry -

1) Resistance :- $I \propto E \therefore E = RI$

(where I = current, E = Potential diff., R = Resistance)
Defⁿ - Current 'I' is directly proportional to E i.e. potential difference betⁿ 2 conductors through which current 'I' is flowing.

Unit :- Ohm . Ω

2) Conductance :- $C = \frac{1}{R}$ Unit $\rightarrow \text{ohm}^{-1} \text{ m}^2$

Defⁿ - It is the reciprocal of resistance (which allows the current to pass through conductors).

3) Specific Resistance or Resistivity :-

$$\rho = R \times \frac{A}{l} \quad (\text{where } \rho = \text{Resistivity}, l = \text{length bet}^n \text{ conductors}, A = \text{area of cross section})$$

Defⁿ - It is the resistance offered by unit length 'l' to the passage of electricity through it.

Unit :- Ohm . cm

4) Specific Conductance or Conductivity :-

$$k = \frac{1}{\rho} = \frac{1}{R} \times \frac{l}{A} = C \times \frac{l}{A} \quad [\text{where } k = \text{specific conductance}, l = \text{length bet}^n \text{ conductors}, A = \text{area of cross sec.}, C = \text{conductance}]$$

Defⁿ - It is the reciprocal of resistivity i.e. the conductance of solⁿ betⁿ 2 parallel electrodes kept 1 cm apart having area of cross sec. 1 cm^2

Unit :- $\text{ohm}^{-1} \cdot \text{cm}^{-1}$ or $\text{mho} \cdot \text{cm}^{-1}$

5) Cell constant :- Unit = cm^{-1}
 $= \frac{1}{A}$ or $= \frac{\text{specific conductance}}{\text{measured conductance}}$

Defⁿ - It is the ratio of specific conductance to the measured conductance.

6) Equivalent Conductance (Λ_v) :-

$$= \frac{k \times 1000}{c} \quad [\text{where } k = \text{sp. conductance}, c = \text{conc. of soln}]$$

Unit :- $\text{ohm}^{-1} \cdot \text{cm}^2 \cdot \text{gm-equiv}^{-1}$

Defⁿ - It is the conductance of a soln containing 1 gm equivalent of an electrolyte at any particular concentration when placed betⁿ two electrodes 1 cm apart.

7) Molar conductance (Λ) :-


$$= \frac{k \times 1000}{M} \quad [\text{where } k = \text{sp. conductance}, M = \text{molarity of soln}]$$

Unit :- $\text{ohm}^{-1} \cdot \text{cm}^2 \cdot \text{gmol}^{-1}$

Defⁿ - It is the conductivity of the soln containing 1 mol of electrolyte when placed betⁿ two large electrodes 1 cm apart.

* Conductivity cell :-

Diagram :-



Construction :-

- i) Measurement of conductance of a solution is done by using conductivity cell with alternating current.
- ii) The most widely used conductivity cell has a pair of Pt electrodes kept 1 cm apart from each other.
- iii) These Pt electrodes/Plates are coated with Pt black which are welded to Pt wires in two respective glass tubes.
- iv) These are connected to Cu wires of circuit dipped in Hg (mercury) filled in glass tubes.
- v) The whole conductivity cell is connected to wheatstone's bridge.

* Conductometric: Acid-Base titrations :-

Types -

- 1) Strong acid vs strong base
- 2) Strong acid vs weak base
- 3) Weak acid vs strong base
- 4) Weak acid vs weak base


1) Strong Acid Vs Strong Base :- (Reⁿ, Titration curve, Explanation)

e.g. HCl vs NaOH

Reaction -



Titration curve -




- Explanation :- i) Initially conductance was observed high because of strong acid (dissociation complete dissociation).
- ii) During titration, as addⁿ of NaOH takes place, conductance decreases rapidly due to replacement of highly mobile H^+ by less mobile Na^+ ions.
- iii) At the end point, conductance is observed due to strong electrolyte NaCl.
- iv) After equivalence/end point, conductance rises due to addⁿ of excess of NaOH (strongly dissociating base),


2) Strong Acid Vs Weak Base →
(Reⁿ, Titration curve, Explanⁿ)



Reaction -



Titration curve -




Explanation -


- i) Initially conductance was observed high because of strongly dissociating HCl.
- ii) On addition of NH₄OH, conductance decreases rapidly due to replacement of highly mobile H[⊕] by slow moving NH₄⁺ ions.
- iii) At the end pt / equivalence point, conductance is observed due to NH₄Cl, electrolyte formation.
- iv) After eq. / end point, conductance remains constant or slightly decreases due to addⁿ of excess of weakly dissociating NH₄OH.

3) Weak Acid vs Strong Base \rightarrow
(Reⁿ, Titration curve, Explanation).
e.g. CH_3COOH vs NaOH

Reaction :-



Titration Curve -



Explanation :-


- i) Initially conductance was observed low because of weakly dissociating CH_3COOH .
- ii) On the addition of NaOH , conductance increases gradually due to formation of sodium acetate ($\text{CH}_3\text{COO}^{\ominus}\text{Na}^{\oplus}$) at an equivalence or end point.
- iii) After equivalence / end point, on excess addⁿ of NaOH conductance rises rapidly due to fast moving OH^{\ominus} & Na^{\oplus} ions.

Page: 2-17
Date: 11


4) Weak Acid Vs Weak Base :- (Reⁿ, Titration Curve, Explanation)

e.g. CH_3COOH vs NH_4OH

Reaction :-



Titration Curve :-



Explanation :-

- Initially conductance was observed low due to weakly dissociating CH_3COOH .
- On addition of NH_4OH from burette, conductance increases slowly due to formation of strong electrolyte ammonium acetate i.e. $\text{CH}_3\text{COO}^- \text{NH}_4^+$ at eq./end point.
- After eq./end point, conductance remains almost constant or slightly decreases due to addition of excess of weakly dissociating NH_4OH soln.

* pH metry :- $pH = -\log [H^+]$ $pH = p =$
Potential of Hydrogen

* Buffer solution -


Defⁿ - The solution which is used to maintain pH of given solution even if small quantity of acid or base is added in it.

Types of Buffer solⁿ - ① Acidic & Basic



Acidic Buffer

Defⁿ It is a mixture of weak acid & salt of weak acid with strong base.




Formula

$$pH \text{ of acidic buffer} = pK_a + \log \frac{[\text{salt}]}{[\text{weak acid}]}$$

Basic Buffer

It is a mixture of weak base & salt of weak base with strong acid.



$$pH \text{ of basic buffer} = pK_b + \log \frac{[\text{salt}]}{[\text{weak base}]}$$


* Standardization or Calibration of pH meter:-

- i) pH meter is standardized / calibrated using pH 4, 7, 9 solⁿ using calomel as a reference & glass as an indicator electrode.
- ii) Each time electrodes are washed with distilled water. Thus, pH meter is calibrated / standardized for acidic, neutral & basic solutions.

* Measurement of pH of Strong Acid vs Strong Base :-


e.g. HCl vs NaOH

Reaction :-



Explanation :-

- Initial stage - Initially burette is filled with std. NaOH solⁿ & pH is noted for given acidic solⁿ.
- Before eq. point - Before addition of NaOH from burette, pH of acidic solution (HCl) is observed as 1 to 2.
- At eq. point - During the addition of NaOH, pH goes on increasing, pH = 7
- After eq. point - After reaction completion, i.e. after equivalence point, pH goes on increasing towards basic value & becomes constant at pH 10.



Calculation - ① Normality of HCl \Rightarrow

$$N_1 V_1 (\text{HCl}) = N_2 V_2 (\text{NaOH})$$

② Strength of HCl \Rightarrow Normality of HCl \times Eq. wt. of HCl

$$1 \text{ ml } 1\text{N NaOH} \equiv 36.5 \text{ mg HCl}$$